

Free Spirit Headset UL-8100 for Clear Vision



Brand: Orgo-Life Quantum MedBed Energy Technology®

Product Code: UL-8100

Availability: 5

Weight: 0.00kg

Dimensions: 0.00cm x 0.00cm x 0.00cm

Price: \$1,280.00

Ex Tax: \$1,280.00

Short Description

Free Spirit Headset UL-8100 for Clear Vision & Vibrational Alignment

Head Meridian Modules with Quantum Push Technology — designed for the Orgo-Life® MedBed Quantum Energy System.

Compatible with: Orgo-Life® Quantum Energy MedBed Technology®

Technology: Quantum Push + Passive PEMF | Non-invasive | 9 vibrational modules (UL-8100) or the corporate 5-module version (UL-8100-A)

Description

PRODUCT DESCRIPTION

Discover the **UL-8100 – Free Spirit Headset**, a mauve vibrational wellness head device designed to support daily well-being and enhance natural clarity of mind. Crafted with soft fabric and integrated velcro straps, it allows the energetic modules to be positioned all around the head, supporting meridian harmony and subtle vibrational alignment.

Each headset includes **5 washable inner head-liners**, perfect for shared use in wellness centers, therapy practices, or families.

The **UL-8100 full kit** contains **9 vibrational modules**, while the **corporate UL-8100-A kit** includes **5 modules**.

Each module has a **positive side and a negative side**, allowing intuitive placement beneath the head or inside the headset depending on the intention of the session.

This is **not a medical device**, but a vibrational wellness tool designed to support mental relaxation, clarity, inner balance, and a natural sense of presence. Every kit includes a **clear and simple instruction guide** to begin immediately.

Ideal for:

- Meditation
 - Relaxation
 - Mental clarity support
 - Energy sessions
 - Wellness centers
 - Home well-being practices

Mission: Reactivate vibrational flow, release old energetic tension, and support the body's natural subtle intelligence.

MISSION OF THE FREE SPIRIT HEADSET

- Support a calmer, more centered inner space
 - Encourage harmonization of the personal energy field
 - Enhance relaxation, meditation, and alignment practices
 - Provide a simple, natural vibrational tool for everyday wellness

Vibrational Headset – Harmony & Clarity

The **Free Spirit Headset UL-8100** is a specialized vibrational tool designed to support the natural energetic balance of the body and promote a calmer, clearer mental presence.

Its mauve design with adaptive velcro allows the vibrational modules to be positioned to encourage meridian flow and energetic alignment.

Each headset includes **5 washable liners**, ideal for personal, family, or professional use.

The UL-8100 comes with **9 vibrational modules**, while the UL-8100-A includes **5**.

These modules have a positive side and a negative side, allowing customized placement beneath the head or within the headset depending on the vibrational

intention.

This is **not a medical device**, but a vibrational wellness tool, fully aligned with **Law 31** standards.

It is compatible with **animals, children, highly sensitive individuals**, and those who desire a gentle, subtle, and professional vibrational approach.

Choosing Orgo-Life® means embracing a modern, natural approach aligned with the laws of living energy—supporting human well-being without replacing human intuition or health professionals.

MAIN VIBRATIONAL BENEFITS

(Non-medical, Law 31 compliant)

- **Meridian harmonization**
Supports smoother vibrational flow around the head.
- **Mental clarity & enhanced presence**
Ideal for refocusing and calming the mind.
- **Deep relaxation**
Encourages a peaceful, stable inner state.
- **Grounding & energetic balance**
Fosters a sense of stability and full-body alignment.
- **Improved inner focus**
Enhances meditation, subtle energy work, and spiritual practices.
- **Daily well-being support**
A simple tool for vibrational renewal throughout the day.
- **Professional or personal use**
Perfect for practitioners, spas, massage therapists, and energy workers.

WHO IS IT FOR?

Recommended for:

- Energy therapists and holistic practitioners
- Vibrational wellness centers
- Individuals who already own an Orgo-Life® MedBed
- Anyone seeking a gentle, non-invasive vibrational approach
- Meditation and breathing practices
- Reiki / energy work
- Conscious breathing sessions
- Vibrational relaxation

- Personal energetic harmonization
- In-office well-being sessions
- Home use for an energetic “reset”

TECHNICAL FEATURES

- Set of 9 vibrational modules using Tesla-inspired micro-coils, designed to emit subtle vibrational fields supporting head and body meridian balance
 - Passive bio-resonance technology — no battery, no current, no electronics
 - Exclusive Orgo-Life® lunar & quantum encoding for a unique vibrational signature
 - Soft, hypoallergenic velcro system for adaptable placement
 - Specialized modules for dense or stagnant vibrational areas
 - Natural synergy with Orgo-Life® Quantum Energy PEMF beds
 - Includes downloadable PDF meridian chart from the Orgo-Life® website

Professional Wellness Use

A fully flexible, reversible, non-invasive vibrational tool — no glue, no drilling, no permanent attachments.

ENHANCED FEATURES

- Soft mauve fabric, breathable, flexible
 - Dual-sided velcro for precise positioning
 - Fully compatible with all FreeSpirit vibrational modules
 - No metal — no interference — pure vibrational clarity
 - Designed for energetic practices: meridian harmonizing, alignment, meditation
 - Universal fit: **40 cm x 40 cm, <90 g**, adapts to 97% of head shapes
 - Includes 5 washable liners (hand-wash recommended)
 - Reusable, lightweight, ideal for personal or professional use

SESSION EXPERIENCE —

As soon as the FreeSpirit modules activate their subtle signature, the body naturally enters a receptive state. Most users describe the first minutes as a gentle vibrational awakening.

You may feel:

- Fine tingles, like tiny sparks beneath the scalp
 - A soft, comforting warmth around the head
 - Pleasant heaviness releasing accumulated tension
 - Subtle emotional loosening, deep breaths, soft exhalations
 - A bubble of stillness forming around you
 - Circular pulsations (“energy rings”) rotating gently
 - Mental realignment — thoughts becoming clearer
 - A strong sense of presence and a grounded inner silence
 - A stable, recharged vitality at the end of the session

Every experience is unique — but consistently described as **deep, grounded, and harmonizing**.

ENERGETIC & TRADITIONAL INSPIRATION

- Aligned with Traditional Chinese Medicine meridian principles
 - Inspired by Ayurvedic energetic alignment
 - Needle-free, gentle vibrational stimulation
 - A blend of ancient symbolism and modern vibrational technology
 - Acts as an amplifier for presence, clarity, and energetic awareness

CARE & DURABILITY

- Clean with a dry or lightly damp cloth (avoid immersion)
 - Hand-built in certified Orgo-Life® workshops
 - 10-year manufacturing warranty
 - Professional users may activate extended PRO warranty with a \$25/year license

VIBRATIONAL TESTIMONIALS (Law 31 compliant)

“From the first session, a soft wave moved through my whole body. My mind felt like someone opened a window inside.”

— **Julie C., Holistic Assistant (Abitibi)**

“I use my mind a lot for work. This headset feels like pressing a mental reset button. I finish grounded, clear, calm.”

— **Michel D., Consultant, Ontario**

“It feels like technology from another timeline... reconnecting me to my own thoughts and my center.”

— **Claire L., Massage Therapist, Sherbrooke**

“I felt circular vibrations around my head — it brought me into meditation instantly.”

— Romain F., Energy Worker, New Brunswick

PRICING & AVAILABILITY

Limited quantities — each module is individually hand-calibrated.

Ships from Québec. International delivery available.

You will receive an email confirming the exact shipping date.

Estimated delay: **7–15 days** before departure from Québec.

TECH SPECS SUMMARY

WEIGHT KIT: approx. **3.3 kg**

MODULES WEIGHT: **175–245 g each**

INCLUDED: 9 module with polarity, mauve headset, 5 liners, laminated chart (EN/FR), PDF instructions

MODEL: **UL-8100 PRO – 9 vibrational modules**

CORPORATE MODEL: **UL-8100-A – 5 modules**

TECHNOLOGY: passive Quantum Push

PACKAGING BOX: **16 x 16 x 5 inches**, approx. **3 kg / 6 lbs**

POLARITY: positive and negative sides on each module

WARRANTY: 10 years with \$25 annual license

BARCODE UL-8100: **787416647455**

BARCODE UL-8100-A: **787416648322**

Note:

Your Orgo-Life® device is fully handmade. Minor bubbles, micro-holes or imperfections may appear when viewed from 5–10 cm — this is normal and part of the artisanal process.

This is not a medical device. For any medical condition, consult a qualified healthcare professional.

METHOD OF USE

(Free Spirit Headset UL-8100 – Well-Being, Law 31 compliant)

Preparation & Alignment – Short Premium Version (before the methods)

To fully benefit from the Free Spirit Headset, a small preparation can transform the session into a deeply soothing and spiritually aligned moment.

The goal is simple: **create a clear, calm, and receptive inner space.**

Recommended preparation (at home or in a wellness center)

- Take a warm bath or a quick shower to release body tension.
- Slightly air out the room to ensure good humidity and a feeling of fresh air.
- Turn off distractions: phone, notifications, noise.
- Sit or lie down in a stable position: lying, seated, or semi-reclined.
- Place a light blanket if you like to feel enveloped.
- Choose a simple intention: clarity, peace, grounding, reset.

Spiritual alignment (short version)

Sit or lie down.

Place your hands on your heart or simply alongside your body.

Take 3 deep breaths.

Imagine your head as an inner lighthouse gently lighting up.

Let the headset amplify this perception and guide your presence toward a clearer, more centered space.

Conscious breathing

- Inhale for 4 seconds
 - Pause for 2 seconds
 - Exhale for 6 seconds
- Repeat 7 to 10 cycles

This breathing synchronizes perfectly with the vibrational modules and helps “empty the mind” — not by forcing, but by letting thoughts drift like clouds.

Mental clarity – mini-ritual

1. Close your eyes
2. Let the headset settle
3. Think of one single phrase:
“I return to my center”
4. Observe your head becoming lighter
5. Let the modules create a clear inner space

For wellness centers

- Install the person with soft music or complete silence
 - Use the washable under-headsets between each client
 - Begin with guided breathing
 - Allow 3 to 5 minutes of integration at the end of the session
 - Always remind that this is a vibrational, non-medical tool (Law 31)

METHODS

1. Mental Clarity Method —

Place the modules around the forehead, temples, and upper skull.

Encourage slow breathing for a few minutes.

Users generally report a sensation of “mental clearing,” like fresh air entering the mind.

Ideal for clearing mental fog, starting a meditation, or preparing for creative focus.

2. Meridian Harmonization Method —

Arrange the modules in a complete circle around the head, forming a “vibrational crown.”

Remain still with eyes closed.

This method helps encourage a more fluid energetic flow and enhances the sensation of overall coherence.

Perfect for grounding, recalibrating, or inner balancing techniques.

3. Deep Relaxation Method —

Place the main module at the top of the head and distribute the others along the back of the skull.

This configuration gently releases accumulated tension and promotes deep unwinding.

Recommended before sleep, after long workdays, or during emotional overload.

4. Guided Meditation Method —

Align the modules along the central axis: *forehead ? top ? back of head.*

This creates a stable vertical energetic line that supports mindfulness and inner stillness.

Ideal for meditation teachers, breathwork, and introspective sessions.

5. Rapid Energetic Reset Method —

Position the modules freely based on intuition or momentary needs.

Sessions of 3–7 minutes can reset internal order, reduce agitation, and restore calm clarity.

Useful for transitions between tasks, grounding after stress, or quick recalibration.

6. Professional Vibrational Care Method —

For practitioners and wellness centers:

Create a vibrational crown around the head, then place one module beneath the head on a soft cushion.

This configuration enhances the receptivity of any energetic, vibrational or holistic care session.

7. Children / Hypersensitivity Method —

Use only 3 to 5 modules for a gentler vibrational flow.

Adapted for children, highly sensitive individuals or anyone new to vibrational tools.

Always supervise children during use.

Provides calmness, grounding, and soft energetic support.

Product Gallery





Orgo-Life® Product Code : **UL-8100**

10 years Limited Warranty with annual operation license
Garantie Limitée de 10 ans avec les droits de licence annuel



Free Spirit Select Wellness Series P80 with Orgo-Life® Quantum MedBed Energy Technology
Série Cosmos Free Spirit et Bio-dôme P80 avec technologie Orgo-Life® MedBed Quantique

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

10 years Limited Warranty with annual operation license
Garantie Limitée de 10 ans avec les droits de licence annuel



Free Spirit Select Wellness Series P80 with Orgo-Life® Quantum MedBed Energy Technology
Série Cosmos Free Spirit et Bio-dôme P80 avec technologie Orgo-Life® MedBed Quantique

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

10 years Limited Warranty with annual operation license
Garantie Limitée de 10 ans avec les droits de licence annuel



Free Spirit Select Wellness Series P80 with Orgo-Life® Quantum MedBed Energy Technology
Série Cosmos Free Spirit et Bio-dôme P80 avec technologie Orgo-Life® MedBed Quantique

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

10 years Limited Warranty with annual operation license
Garantie Limitée de 10 ans avec les droits de licence annuel

Included:
- 4 Undercap Ring
- 1 sensor Head cap
- 8 head modules

Inclus:
- 4 Doubles anneaux-organes
- 1 capteur de tête module
- 8 modules de tête



Free Spirit Select Wellness Series P80 with Orgo-Life® Quantum MedBed Energy Technology
Série Cosmos Free Spirit et Bio-dôme P80 avec technologie Orgo-Life® MedBed Quantique

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

10 years Limited Warranty with annual operation license
Garantie Limitée de 10 ans avec les droits de licence annuel



Free Spirit Select Wellness Series P80 with Orgo-Life® Quantum MedBed Energy Technology
Série Cosmos Free Spirit et Bio-dôme P80 avec technologie Orgo-Life® MedBed Quantique

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

10 years Limited Warranty with annual operation license
Garantie Limitée de 10 ans avec les droits de licence annuel



Free Spirit Select Wellness Series P80 with Orgo-Life® Quantum MedBed Energy Technology
Série Cosmos Free Spirit et Bio-dôme P80 avec technologie Orgo-Life® MedBed Quantique

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

Orgo-Life® MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Orgo-Life® MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Orgo-Life® MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology



Free Spirit MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Free Spirit MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Free Spirit MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology

MedBedEnergy.com Orgo-Life.com Med-Bed.ca

Orgo-Life® Product Code : **UL-8100**

Orgo-Life® MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Orgo-Life® MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Orgo-Life® MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology



Free Spirit MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Free Spirit MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology
 Free Spirit MedBed Series PBD with Orgo-Life® Gasoline MedBed Energy Technology

MedBedEnergy.com Orgo-Life.com Med-Bed.ca